The natural way

SISSEL® NECK PILLOWS

- | Gebrauchsanweisung
- | Instructions for use
- | Mode d'emploi
- | Gebruiksaanwijzing
- | Instrucciones de uso
- | Istruzioni per l'uso

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SISSEL® NECK PILLOWS

General information

Read the following user manual carefully and observe the safety instructions and warnings. The SISSEL neck pillow is an orthopedic pillow that supports the cervical spine due to its special shape and thus relieves the entire shoulder-neck area. The ergonomically shaped neck pillow SISSEL Compact Pillow combines the functionality of the proven SISSEL neck pillows with compact dimensions suitable for travelling. So you can rely on your usual sleeping comfort even when you are on the road. The product and accessories may only be used in accordance with these operating instructions or as instructed by appropriately trained specialists. SISSEL GmbH accepts no liability for the use of the system outside the applications described in these instructions for use. The user manual must be kept carefully and within easy reach, and must be enclosed if the product is passed on. The regulatory requirements applicable to this product in the respective country of use must be taken into account.

Safety instructions

- ! Check the undamaged condition of the product before each use. In case of damage, the function of the product and safety is no longer guaranteed.
- ! Do not use the product in the close area of heat sources (e.g. radiators, stove top, open flame, flying sparks) and avoid direct sunlight.
- ! Keep the area of use free of sharp or pointy objects.
- ! Use only accessories recommended by the manufacturer.
- ! In case of uncertainty or existing health problems in the shoulder neck area or cervical spine, use the product only after consulting a doctor.

Product information

Scope of delivery

consisting of one of the following products:

- 1 x SISSEL Soft Bambini, 1 x instructions for use, 1 x cover
- 1 x SISSEL Soft, 1 x instructions for use, 1 x cover
- 1 x SISSEL Soft Plus, 1 x instructions for use, 1 x cover
- 1 x SISSEL Soft Deluxe, 1 x instructions for use, 1 x cover
- 1 x SISSEL Classic, 1 x instructions for use, 1 x cover
- 1 x SISSEL Classic Plus, 1 x instructions for use, 1 x cover
- 1 x SISSEL Compact Pillow, 1 x instructions for use, 1 x cover

Technical Data

Dimensions (L x W x H):	SISSEL Soft Bambini:	approx. 37 x 25 x 9 cm
	SISSEL Soft:	M: approx. 47 x 33 x 11 cm L: approx. 47 x 33 x 14 cm
	SISSEL Soft Plus:	approx. 47 x 33 x 11 / 14 cm
	SISSEL Soft Deluxe:	approx. 63 x 33 x 10 / 12 cm
	SISSEL Classic:	M: approx. 47 x 33 x 11 cm L: approx. 47 x 33 x 14 cm
	SISSEL Classic Plus:	approx. 47 x 33 x 11 / 14 cm
	SISSEL Compact Pillow:	approx. 37 x 25 x 9 cm

Material:	SISSEL Soft Bambini:	Filling: Viscoelastic foam Cover: 100% cotton
	SISSEL Soft:	Filling: viscoelastic foam Cover: Velour (75% cotton, 25% polyester)
	SISSEL Soft Plus:	Filling: Viscoelastic foam Cover: Velour (75% cotton, 25% polyester)
	SISSEL Soft Deluxe:	Filling: Viscoelastic foam Cover: Satin (50% cotton, 50% polyester)
	SISSEL Classic:	Filling: Polyurethane foam Cover: Velour (75% cotton, 25% polyester)
	SISSEL Classic Plus:	Filling: Polyurethane foam Cover: Velour (75% cotton, 25% polyester)
	SISSEL Compact Pillow:	Filling: Viscoelastic foam Cover: Velour (75% cotton, 25% polyester)
Weight:	SISSEL Soft Bambini:	approx. 0,60 kg
	SISSEL Soft:	M: approx. 0,98 kg L: approx. 1,16 kg
	SISSEL Soft Plus:	approx. 1,16 kg
	SISSEL Soft Deluxe:	approx. 1,20 kg
	SISSEL Classic:	M: approx. 0,58 kg L: approx. 0,74 kg
	SISSEL Classic Plus:	approx. 0,74 kg
	SISSEL Compact Pillow:	approx. 0,60 kg

Note: The SISSEL Soft Plus, the SISSEL Soft Deluxe and the SISSEL Classic Plus are height-variable by an insert to adjust the neck pillow individually optimal.

Intended use

Use the SISSEL neck pillow instead of your usual pillow for sleeping. The pillow is designed especially for side sleepers, but also for back sleepers. The special shape positions the cervical spine in an anatomically correct, physiological position and relieves the neck muscles and intervertebral discs. In this way, the SISSEL neck pillow makes a significant contribution to relaxing the neck muscles and thus promoting deep, healthy, restful sleep. Lie on your side with your head on the SISSEL neck pillow so that the high side is against your shoulder and the cervical spine is supported by the pad. The SISSEL neck pillow is correctly positioned when the high side fills the area where the cervical spine meets the shoulder. When using the SISSEL neck pillow, the new sleeping position may initially cause a change in sleeping habits for a few nights. As a rule, this adjustment lasts a few days. Only in exceptional cases can it take longer.

Indikation

- Pain in the shoulder-neck area and the cervical spine
- Tension in the neck muscles with accompanying tension headaches
- Arthrosis and / or intervertebral disc problems of the cervical spine (consultation with specialist)
- Therapeutic support for cervical spondylosis (consultation with specialist)
- Accompanying therapy for chronic instability syndrome of the cervical spine (consultation with specialist)
- Rehabilitation support for acute or chronic neck pain, tension of the neck muscles and shoulder girdle (consultation with specialist)
- Prevention of tension and overstraining malpositions through anatomically correct positioning of the cervical spine during sleep
- Avoidance of overloading the cervical spine through physiological positioning of the head in the side position during sleep

Contra Indication

- Lack of compliance (due to fear or pain)
- Acute injuries and / or sprains of the cervical spine
- Acute cervical spine syndrome (cervical syndrome)
- Acute cervical vertigo (dizziness)
- Acute radiating pain
- Discomfort
- Movement disorders
- Acute stroke
- Acute herniated disc of the cervical spine
- Acute inflammatory processes such as meningitis and encephalitis
- Acute or chronic nerve root inflammation of the cervical spine
- Spastic torticollis (spasmodic torticollis)
- Neurological diseases or malformations
- Tumour diseases

Regulatory requirements

The current product of SISSEL GmbH complies with the following regulatory requirements:

- Medical Devices Regulation (EU) 2017 / 745

 All serious adverse events that occur in relation to the product must be reported immediately to the manufacturer and the competent authority.

Accessories

The product may only be used with the accessories supplied.

Usage

1. Check the packaging and the product for any damage upon delivery.

 Never use a damaged product.

2. Check the product for integrity and contamination before each use.

3. Use the product as intended.

Note: Cleaning measures may be necessary after use (see Care & hygiene).

Care & hygiene

This product is suitable for repeated use. In general, the applicable national regulatory requirements for the hygienic reprocessing of medical devices must be observed and fault-free function must be ensured. SISSEL neck pillows cannot be disinfected. If used by more than one person, additional hygiene measures may be necessary. The PU foam core and the core made of visco-elastic foam can be hand-washed at lukewarm temperature using a mild detergent. Carefully squeeze out by hand, never wring out! Then leave to air dry on a grid (place a towel underneath) for a sufficiently long time (approx. 1 week). Do not dry in the sun or on the heater. Discoloration of the pillow core is possible and does not affect the function of the product. The velour cover and the cotton cover are machine washable at 40 °C. The satin cover can also be washed at 60 °C. Then air dry. Follow the washing label.

Storage: For storage, protect the product from direct sunlight and store it in a clean, dark and dry place when not in use.

Disposal & environment

The product and the accessories supplied with it can be disposed of with household waste. No special types of packaging are used. The respective valid disposal guidelines must be obeyed.

Warranty

The warranty is governed by the legal regulations. Any warranty is excluded if goods delivered by SISSEL GmbH have been processed, treated or modified by others without consent or if the instructions for use have not been followed. The user of this product must comply with the regulatory requirements for medical devices applicable at the place of use.

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